

Beamline 5.4 is presently under construction.

It will be located on top of the LINAC/BTF cave, inside the main ALS storage ring!

Source Bend magnet synchrotron radiation

Acceptance Angle 17 mrad (vert) x 69 mrad (horiz)

Scientific Applications Biological samples, environmental samples, novel compounds, forensic stud

Endstations

Beamline 5.4.2 Beamline 5.4.3

Beamline 5.4.1

Description□ Infrared Spectromicroscopy to be determined High resolution Far-IR to mid-IR spectroscopy

-1

Frequency Range 600-10,000 cm

-1 20-5,000 cm

Spot Size
3-10 um (diffraction limited)
1 mm

Resolution

 $0.125 \ \text{cm}$ $^{-1}$ $0.00096 \ \text{cm}$ $^{-1}$

